

# Assessing Antibody Effects on Platelet Function: *In Vitro* Assays for De-Risking in Antibody Discovery and Development




C. Jewell, W. Tomlinson, A. Johnson, P. Cato,  
S. Bhaker, N. Dovlatova. Platelet Services, Nottingham, UK

## Background


- Therapeutic antibodies (mAbs, ADCs, bispecifics, fragments) dominate new drug pipelines. Some antibodies and IgG-containing immune complexes can activate platelets via FcγR1a, increasing thrombotic risk.
- We evaluated which *in vitro* platelet assays are most applicable for detecting antibody- and immune complex-mediated effects on human platelet function.**

## Mechanisms of Antibody Effects on Platelets





Human platelets express **FcγR1a**. When an antigen is present, **IgG** antibodies such as anti-CD9, bind and **crosslink FcγR1a** leading to platelet activation.

**IgG-containing immune complexes**, such as anti-CD40 with soluble CD40L, can also activate platelets through FcγR1a signalling.



## Methods

- Blood was collected from healthy human volunteers.
- Platelet activation was measured in whole blood (WB) and washed platelets (WP) by flow cytometry (P-selectin expression).
- Platelet aggregation was measured in platelet-rich plasma (PRP) and WP by changes in light absorbance over time.
- All assays were carried out on 96-well plates and measurements were performed in duplicate.



1

## Results: Effects of Anti-CD9, Anti-CD226, and Anti-CD47 on Aggregation in Platelet Rich Plasma




Figure 1. Aggregation assessed in PRP with anti-CD9, anti-CD226 and anti-CD47 (n=7-9; mean ± SEM).

- Volunteers were classified as **high** or **low** responders based on aggregation with 1 µg/mL of anti-CD226: >80% aggregation = "**high responders**" (Based on data from 36 volunteers).


- The three antibodies induced aggregation with varying potency; relative potency was consistent between **high** and **low** responders.

2


## Results: Anti-CD226-Induced Platelet Responses Across Platelet Assays and Matrices

### Platelet Activation

#### Whole Blood




#### Washed Platelets



### Platelet Aggregation

#### Platelet-Rich Plasma



#### Washed Platelets




Figure 2. Activation (WB, WP) and aggregation (PRP, WP) co-stimulated with ADP and anti-CD226 in **high responders** (n=5-6; mean ± SEM).

3

## Results: Anti-CD40L:sCD40L Immune Complexes Induce Platelet Aggregation




Figure 3. Anti-CD40L:sCD40L immune complex effects on platelet aggregation in WP obtained from **high responders** (n=5-7; mean ± SEM).

- Ruplizumab (an anti-CD40L), discontinued in clinical development due to thromboembolic risk (Kawai *et al.*, 2002), forms immune complexes (ICs) with soluble CD40L (sCD40L).

- These ICs induced concentration-dependent platelet aggregation in washed platelets, with the highest response observed at a 1:3 molar ratio. In contrast, no clear aggregation was observed in platelet-rich plasma (data not shown).

## Conclusion

- Platelet activation and aggregation **assays performed in a 96-well plate format** can clearly demonstrate antibody-induced platelet responses; **anti-CD226** showed **concentration-dependent effects**, most pronounced in **washed platelets**.
- Platelet reactivity via **FcγR1a** varies substantially among human volunteers; **anti-CD9, anti-CD226 and anti-CD47** showed **volunteer-dependent responses**, suggesting that recruiting "**high responders**" improves the detection of secondary pharmacology effects of antibodies on platelets.
- Ruplizumab (anti-CD40L) and soluble CD40L IC**, induced **platelet aggregation**, reflecting the proposed FcγR1a-mediated mechanism underlying its reported thromboembolic risk.
- In vitro* platelet assays can support early de-risking during antibody discovery and development.**

info@plateletservices.com

